

# AMBIC®

## EasiDoser™



### OPERATING INSTRUCTIONS

MODE D'EMPLOI



ANLEITUNG

*leading best practice in livestock health management*

## Contents

Introduction .....	3
Safety .....	3
Specifications .....	3
Installation (refer to Fig.1, 2, 4) .....	4
Layout .....	4
Electrical Connections .....	5
Opto Isolator Connection .....	6
Initial Set Up – Electrical Connections & Calibration/Programming Functions .....	7
EasiDoser™ – Power Up/Down & Normal Run Displays .....	7
EasiDoser™ – Basic Menu Option Displays .....	8
EasiDoser™ – Set Two Pump Operations .....	8
EasiDoser™ – Changing Pin Code .....	10
EasiDoser™ – Pump 1 Calibration Method .....	12
EasiDoser™ – Pump 2 Calibration Method .....	12
EasiDoser™ – Pump Statistics .....	13
EasiDoser™ – Warning Displays .....	13
Initial Set Up – Priming & Setting up the system ready to use .....	14
Operation .....	14
Maintenance .....	14

# Ambic Equipment Ltd – EasiDoser™ Operating Notes

## Introduction

The EasiDoser™ is an advanced microprocessor-controlled dosing system for dispensing a precise amount of chemical. The dispenser features an electronic circuit board with an LCD display. Programming is set using the keypad on the front of the unit and has a special pump 'Lock-out' that will prevent consecutive dispensing. The EasiDoser™ also has a "delay pump start" feature and a direct mode (Relay Mode) which allows for a timed signal from a microprocessor to control the pump run timing.

The EasiDoser™ is offered in various configurations: one, two or three low volume pumps, or an industrial pump version.

Depending how the unit is programmed, the EasiDoser™ will activate when the circuit receives a 12-230V AC or 12-24V DC machine signal through the Opto Isolator, or closed contact to one of the trigger inputs. Also, it can be used for manual dosing and this is activated by pressing a button on the keypad.

## Safety

The EasiDoser™ Chemical Dosing System is designed exclusively for dosing chemicals by function of Time, Time and Speed or Direct (Relay) mode in unrestricted pressure tanks. Any application outside the use described in this operating manual will be taken to be not in accordance with the intended purpose. The manufacturer/supplier will not be held responsible for any losses arising as a result of such use. The user will take full responsibility for use.

USE IN ACCORDANCE WITH THE INTENDED PURPOSE ALSO INCLUDES COMPLYING WITH THE OPERATING MANUAL AND THE CONDITIONS FOR INSPECTION AND MAINTENANCE.

## Specifications

Power Source	230V – Pumping unit – Fuses	<ul style="list-style-type: none"><li>• 220-240V AC 50Hz</li><li>• T1.0A</li></ul>
	115V – Pumping unit – Fuses	<ul style="list-style-type: none"><li>• 110-120V AC 60Hz</li><li>• T2.0A</li></ul>
Machine Signal -		<ul style="list-style-type: none"><li>• 12 – 230V AC</li><li>• 12 – 24V DC</li><li>• Contact switch</li></ul>
Maximum Length of Distribution Tubing		<ul style="list-style-type: none"><li>• 10 Metres</li></ul>
Operating Temperature		<ul style="list-style-type: none"><li>• 5 – 40 °C</li></ul>
Float Switch		<ul style="list-style-type: none"><li>• Normally Open type</li></ul>

## Installation (refer to Fig.1, 2, 4)

**Electrical Installation** - should ideally be carried out by a Qualified Electrician – The 3-core power supply cord **MUST** remain accessible when the unit has been installed and **MUST** be connected to Earth. THIS ELECTRICAL SUPPLY SHOULD BE PROTECTED BY AN APPROPRIATE EARTH LEAKAGE CIRCUIT BREAKER.



**CAUTION:** To avoid severe or fatal shock, always disconnect main power when servicing the unit.

**Chemical** – is sucked up from a chemical container by peristaltic pump(s). The enclosure housing the EasiDoser™ pump(s) is rated at IP44 (IEC 60529); the power unit (**A**) must, therefore, be positioned in a dry location and close to a 220-240V AC power socket (Maximum Power consumed: 50W).

The **enclosure (A1)** should be positioned no more than 2 metres vertically above the floor and no more than 3 metres distant from the chemical container(s). Use the drilling template to mark out screw positions and fix the enclosure to a suitable flat surface, using the screws and plastic plugs supplied.

The **chemical container(s) (D)** should be placed securely on the floor and the inlet tube(s) (**C**) from the pump enclosure should be connected to the inlet nipple(s) of the peristaltic pump(s) – each PVC inlet tube should be measured for correct length (cutting the excess tube) and warmed before being pushed over the peristaltic pump spigot and secured with the stainless steel hose clips. Make sure that clips are fitted so that the barb of the fitting on the pump is central in the clip and do NOT over tighten the clip (you risk breaking the plastic nipple). Drop the filter end of the inlet tube(s) into the chemical container(s) – ALWAYS ENSURE that each inlet tube has a FILTER FITTED.



**CAUTION:** Wear protective clothing and eyewear when dispensing chemicals or any other materials. Observe safety handling instructions (MSDS) of chemical manufacturers.

## Layout




Fig. 1 - General Layout for Twin EasiDoser™




Fig. 2 – General layout for a Single EasiDoser™

## Electrical Connections

**WARNING!**


**DISCONNECT UNIT COMPLETELY FROM MAINS POWER BEFORE OPENING**

Unscrew the 2 screws in the corners of the enclosure. Hinge lid open to reveal the PCB inside lid.

Terminal blocks may be pulled off from the PCB to allow easy connection of wires into Screw Terminals.

Positions of wires are noted **LEFT >> RIGHT** (see PCB Layout diagram below right)

Figure 3 – EasiDoser™ Switch Operations



**MAINS** = 220-240V AC

- L** = LIVE – **Brown** wire
- E** = EARTH – **Green/Yellow** wire
- N** = NEUTRAL – **Blue** wire

- P1** = Pump 1 – **Pink** wire + **Grey** wire -
- P2** = Pump 2 – **Pink** wire + **Grey** wire -
- P3** = Pump 3 – **Pink** wire + **Grey** wire -


**Opto 4** = Signal wire Pump 1  
 12-230V AC via Opto Isolator  
 12-24V DC via Opto Isolator  
 Close contact direct input

**Opto 5** = Signal wire Pump 2  
 12-230V AC via Opto Isolator  
 12-24V DC via Opto Isolator  
 Close contact direct input

**Opto 3** = Signal wire Pump 3  
 12-230V AC via Opto Isolator  
 12-24V DC via Opto Isolator  
 Close contact direct input


**FS1 & 2** = Float Switch 1 & 2  
 – **Brown** wire  
 – **White** wire

Figure 4 – EasiDoser™ Electrical Connections



## Opto Isolator Connection

For volts-free signal (switch) there is no requirement for Opto Isolator and the switch can be connected directly to the inputs marked on the PCB as Opto 4 & Opto 5. For a 12-24V DC signal or 24-230V AC signal unit, the signal **MUST** be connected via the Opto Isolator(s) supplied with the unit.



**Figure 5 – Opto Isolator**

**Opto Isolator 1**


Brown & Blue – connect to 12-24V DC or 24-230V AC signal  
Red & Black – connect to Opto 4

**Opto Isolator 2**

Brown & Blue – connect to 12-24V DC or 24-230V AC signal  
Red & Black – connect to Opto 5

**Opto Isolator 3**

Brown & Blue – connect to 12-24V DC or 24-230V AC signal  
Red & Black – connect to Opto 3




**Figure 6 – Opto Isolator position on PCB**

The Red & Black wires should be connected to Opto 4 or Opto 5.

Opto 4 triggers the left hand pump and Opto 5 triggers the right hand pump on a 2 pump system.

On a 3 pump system, Opto 3 triggers pump 3.



## Initial Set Up – Electrical Connections & Calibration/Programming Functions

The EasiDoser™ control box is connected up as shown in Figure. 4 on page 5. It includes many useful safety and programmable features, for which the following instructions should be read through carefully to ensure efficient operation of the unit.

### EasiDoser™ – Power Up/Down & Normal Run Displays



































































Action	Display	Comment /Action
Switch Mains Power ON Press on/off Button	<p>Display Number 1</p> <div style="border: 1px solid black; padding: 5px; text-align: center;">                     Ambic Equipment EasiDoser V1.1                 </div>	Displays for 5 seconds. Automatically changes to NEXT display.
	<p>Display Number 2</p> <div style="border: 1px solid black; padding: 5px; text-align: center;">                     Waiting for Trigger ←0→                 </div>	Pump(s) NOT running. No trigger signal applied or start dose button pressed.
Press On/Off Button (when switching off the unit/ leaving on standby)	<p>Display Number 3</p> <div style="border: 1px solid black; padding: 5px; text-align: center;">                     EasiDoser Power Down                 </div>	Displays for 2 seconds, then screen blank. Unit in “Standby / Off” mode.
Press On/Off Button	<p>Display Number 1</p> <div style="border: 1px solid black; padding: 5px; text-align: center;">                     Ambic Equipment EasiDoser V1.1                 </div>	Displays for 5 seconds. Automatically changes to display 2.
Disconnect from Mains Power	<p>Display Number 4</p> <div style="border: 1px solid black; padding: 5px; text-align: center;">                     Powering Down Bye                 </div>	Unit completely “Off” Warning “Bleep” before screen blank.

The EasiDoser™ is supplied as either a single, twin or three pump system and can be set to run up to the maximum number of pumps fitted. This is controlled by the user-programmable settings selected from the menus detailed below.

Pumps are driven by low voltage DC motors and their time or speed of operation is controlled by the voltage applied to each pump motor in order to ensure that the pumps are dispensing the correct quantity of chemical(s).

## EasiDoser™ – Basic Menu Option Displays










### EasiDoser™ – Set Two Pump Operations

Action	Display	Comment /Action
Press  Button twice	<p>Display Number 5</p> <div style="border: 1px solid black; padding: 5px;">                     Enter PIN Code                      —                 </div>	<b>Enter the current PIN Code</b> by pressing the buttons (in sequence):- (default factory PIN =     )
Correct PIN code entered	<p>Display Number 6</p> <div style="border: 1px solid black; padding: 5px;">                     Number of Pumps                      = 2                         </div>	Press   buttons to change the number of pumps fitted.
<b>From Display 6</b> Press  Button once	<p>Display Number 7</p> <div style="border: 1px solid black; padding: 5px;">                     Pump Types                      12VOLT                          </div>	Press   buttons until correct voltage & pump type displayed.
<b>From Display 7</b> Press  Button once	<p>Display Number 8</p> <div style="border: 1px solid black; padding: 5px;">                     Pump 1 Mode                      Manual Mode                          </div>	Press   buttons until correct pump mode for Pump 1 is displayed. (Manual, Signal or Direct)
<b>From Display 8</b> Press  Button once	<p>Display Number 9</p> <div style="border: 1px solid black; padding: 5px;">                     Pump 1 Duration                      00:15                          </div>	Press   buttons until correct pump duration (run time) is displayed. Run time can be set from 00:00 seconds to 05:00 minutes. Default setting: 00:15
<b>From Display 9</b> Press  Button once	<p>Display Number 10</p> <div style="border: 1px solid black; padding: 5px;">                     Pump 1 Start                      Delay = 5s                          </div>	Press   buttons until correct pump Start Delay is displayed. Pump start delay time can be set from 1 second to 60 seconds. Default setting: 5s
<b>From Display 10</b> Press  Button once	<p>Display Number 11</p> <div style="border: 1px solid black; padding: 5px;">                     Pump 1 Inhibit                      00:00                          </div>	Press   buttons until correct pump duration (run time) is displayed. Pump inhibit time can be set from: 00:00 seconds to 30:00 minutes. Default setting: 00:00
<b>From Display 11</b> Press  Button once	<p>Display Number 12</p> <div style="border: 1px solid black; padding: 5px;">                     Pump 1 Speed                      = 80%                          </div>	Press   buttons until correct pump Start Delay is displayed. Pump speed can be set from 40% to 100% Default setting: 80%
<b>From Display 12</b> Press  Button once	<p>Display Number 13</p> <div style="border: 1px solid black; padding: 5px;">                     Conductivity                      Pump1 OFF                          </div>	Press   buttons until correct conductivity is set. Available settings are from 0.5mS to 9mS. Default setting: OFF
<b>From Display 13</b> Press  Button once	<p>Display Number 14</p> <div style="border: 1px solid black; padding: 5px;">                     Pump 1 Float                      N/Open                          </div>	Press   buttons until correct Float Switch is displayed. Default setting: N/Open





**EasiDoser™ – Changing Pin Code**

Action	Display	Comment /Action
Use MENU to select	<p style="text-align: center;">Display Number 21</p> <div style="border: 1px solid black; padding: 5px; text-align: center;">                     Change PIN Code ↑                      0 ↓                 </div>	To change the PIN code used to access certain MENU items.
Press  Button once	<p style="text-align: center;">Display Number 5</p> <div style="border: 1px solid black; padding: 5px; text-align: center;">                     Enter PIN code                      -                 </div>	<b>Enter the current PIN Code</b> by pressing the buttons (in sequence):- (default factory PIN =     )
If incorrect PIN entered	<p style="text-align: center;">Display Number 26</p> <div style="border: 1px solid black; padding: 5px; text-align: center;">                     Incorrect PIN                      Code entered                 </div>	Reverts to <b>Display 21</b> after a short delay,
After correct PIN is entered and a short delay	<p style="text-align: center;">Display Number 27</p> <div style="border: 1px solid black; padding: 5px; text-align: center;">                     Enter new PIN                      ↑ → ↑ → _                 </div>	<b>Enter a new 4-digit PIN Code</b> by pressing any combination of these four buttons:-     New PIN is displayed as it is entered.
After final PIN digit entered and a short delay	<p style="text-align: center;">Display Number 28</p> <div style="border: 1px solid black; padding: 5px; text-align: center;">                     Saved PIN Code                      ↑ → ↑ → _                 </div>	Displays new PIN. Reverts to <b>Display 21</b> after a short delay,

## EASIDOSER™ – CHANGING MENU LANGUAGE

By using the MENU screens, it is possible to set the screen language to either

ENGLISH, FRANCAIS, or DEUTSCH as may be most appropriate.

Use MENU to select

Display Number 5

Enter PIN code

—

The default PIN code is entered by pressing the buttons (in sequence):-

↑ ← → ↓

### SETTING ENGLISH MENU

**From Display 21**

Press ↓ Button once

Display Number 22

Select Language

English ← 0 →

Press → button once to move to SET FRENCH MENU - **Display LF** (see below).

**Display LS**

Press □ Button once confirms English

EasiDoser ↑

Change Language ↓

Press ↓ button once to return to OPERATION MENU (**Display 7 or 8**)

**ENGLISH MENU** set .

**Display E**

### SETTING FRENCH MENU

**From Display LS**

Press → Button once

Select. Langue

Francais ← 0 →

Press → button once to move to **Display LD** or button ← to return to **Display LS**.

**Display LF**

Press □ Button once

EasiDoser ↑

Select. Langue ↓

Press ↓ button once to return to OPERATION MENU (**Display 7 or 8**)

**FRENCH MENU** set.

**Display F**

### SETTING GERMAN MENU

**From Display LS**

Press → Button TWICE

Wahlen Sprache

Deutsch ← 0 →

Press ← button once to return to **Display LF** button ← TWICE to return to **Display LS**.

**Display LD**

Press □ Button once

EasiDoser ↑

Wahlen Sprache ↓












Press ↓ button once to return to OPERATION MENU (**Display 7 or 8**)

**GERMAN MENU** set.

**Display D**












## Ambic Equipment Ltd – EasiDoser™ Operating Notes

### EasiDoser™ – Pump 1 Calibration Method



Action	Display	Comment/Action
	Display Number 23	
Press  Button once	Calibrate <span style="float: right;">↑</span> Pump 1            0 <span style="float: right;">↓</span>	To check the dispensing volume of Pump 1 take the calibration vessel and place distribution /delivery tube into the vessel.
	Display Number 5	
Press  Button once	Enter PIN Code -	The default PIN code is entered by pressing the buttons (in sequence):-    
	Display Number 29	
Once correct PIN entered	Press  to Start Calibration	Press button  once - pump 1 runs for 1 minute and then stops. Then the display reverts to next display.
	Display Number 30	
When Pump stops	Volume Dispensed 1 ml            ← 0 →	Check the volume dispensed and using  enter the correct volume dispensed. Then press button  once.
	Display Number 23	
	Calibrate <span style="float: right;">↑</span> Pump 1            0 <span style="float: right;">↓</span>	Press button  once to select Pump 2 Calibration.

For DUAL Pump Operation BOTH Pumps MUST be calibrated.





### EasiDoser™ – Pump 2 Calibration Method

Action	Display	Comment/Action
	Display Number 24	
Press  Button once	Calibrate <span style="float: right;">↑</span> Pump 2            0 <span style="float: right;">↓</span>	To check the dispensing volume of Pump 2 take the calibration vessel and place distribution /delivery tube into the vessel.
	Display Number 5	
Press  Button once	Enter PIN Code -	The default PIN code is entered by pressing the following buttons (in sequence):-    
	Display Number 29	
Once correct PIN entered	Press  to Start Calibration	Press button  once - pump 2 runs for 1 minute and then stops. Then the display reverts to next display.
	Display Number 30	
When Pump stops	Volume Dispensed 1 ml            ← 0 →	Check the volume dispensed and using  enter the correct volume dispensed. Then press button  once.
	Display Number 24	
	Calibrate <span style="float: right;">↑</span> Pump 2            0 <span style="float: right;">↓</span>	Press button  once to save the configuration.

## Ambic Equipment Ltd – EasiDoser™ Operating Notes

To check the Pump Statistics, switch OFF the unit by pressing the button  on the front of the unit. Then press both the down arrow and the  button at the same time.

### EasiDoser™ – Pump Statistics

Action	Display	Comment /Action
Press  Button and Power on/off button	Display Number 31 <div style="border: 1px solid black; padding: 5px; display: inline-block;">                         Pump                    ↑                          Statistics                ↓                     </div>	Press the  button
Press  Button once	Display Number 32 <div style="border: 1px solid black; padding: 5px; display: inline-block;">                         Pump 1: 0.35L        ↑                          0000:03:08        0 ↓                     </div>	Shows amount delivered by Pump 1 on 1 <sup>st</sup> line and on the 2 <sup>nd</sup> line shows the length of time the pump has been running for.
Press  Button once	Display Number 33 <div style="border: 1px solid black; padding: 5px; display: inline-block;">                         Pump 2: 0.35L        ↑                          0000:03:08        0 ↓                     </div>	Shows amount delivered by Pump 2 on 1 <sup>st</sup> line and on the 2 <sup>nd</sup> line shows the length of time the pump has been running for.

### EasiDoser™ – Warning Displays

#### IF LEVEL PROBES ARE FITTED IN CHEMICAL CONTAINERS

When Drum 1 becomes LOW in chemical	Display Number 34 <div style="border: 1px solid black; padding: 5px; display: inline-block;">                         Float Alarm 1                     </div>	Sounder activates and display shown below <b>automatically alternates</b> with normal operational displays.  <b>** Refill/replace container 1 as soon as possible **</b>
When Drum 2 becomes LOW in chemical	Display Number 35 <div style="border: 1px solid black; padding: 5px; display: inline-block;">                         Float Alarm 2                     </div>	Sounder activates and display shown below <b>automatically alternates</b> with normal operational displays.  <b>** Refill/replace container 2 as soon as possible **</b>
If BOTH Drums 1 & 2 become LOW in chemical	Display Number 36 <div style="border: 1px solid black; padding: 5px; display: inline-block;">                         Float 1 &amp; 2 Alarm                     </div>	Sounder activates and display shown below <b>automatically alternates</b> with normal operational displays (e.g. <b>Display 2 &amp; 3 or 4</b> ).  <b>** Refill/replace BOTH containers as soon as possible **</b>

## Ambic Equipment Ltd – EasiDoser™ Operating Notes

If Optional Level Probes (e.g. AFF/200-50-10 – Float Switch Assembly for 20-25 Litre Container) are fitted, LOW level of Chemical activates Float Switch as below:-

- Sounder will operate 3 times every 2-3 seconds.
- Sounder sequence will continue to operate at 1 minute intervals.
- Warning Displays (as shown above) will alternate with the normal run displays.
- When the container is refilled/replaced, the Sounder and Warning Displays stop automatically only when the pump makes its first Dispensing cycle after refill.
- Sounder CANNOT be muted until level of Chemical rises above “LOW”.

### Initial Set Up – Priming & Setting up the system ready to use

Switch on the power at the electric socket and on the front of the pumping unit. Set the program for direct mode. Press the appropriate button in the front of the unit to activate the required pump and liquid will start to be drawn up through the pump(s). Allow unit to run until liquid has passed through the pump(s) and is flowing out in a continuous stream from the end of the distribution pipe. Reset the direct mode to your required setting.

### Operation

Depending on the set up of the unit the pump(s) will either run for a set duration or continuously until the signal is removed:

Manual mode – pressing the dose button for 1 second activates for a set duration.

Signal mode – receipt of either a 12-24V DC or a 12-230V AC signal activates for a set duration.

Direct mode – receipt of either a 12-24V DC or a 12-230V AC signal activates the unit and it will run continuously until this signal is removed.

**WARNING – Refill chemical container(s) BEFORE IT IS EMPTY and SWITCH OFF PUMP UNIT whilst refilling.** In the event that chemical supply is exhausted, it will be necessary to re-prime system as described in the “Initial Set Up section” above.

### Maintenance

In order to preserve the efficiency and reliability of the EasiDoser™ system, we recommend that the following parts be inspected regularly for physical wear and replaced EVERY 1 – 2 YEARS, depending on the chemical used and the frequency of use.

Peristaltic pump tubes - Grease the tubes with silicone every six months

Inlet tube filters (part ATS/419). - Clean regularly with warm water

For optimum results and safety use all the parts supplied by the manufacturer.





Ambic Equipment Limited,

1 Parkside, Avenue Two, Station Lane,  
Witney, Oxfordshire, OX28 4YF. England  
Tel: +44 (0)1993 776555 Fax: +44 (0)1993 779039

[www.ambic.co.uk](http://www.ambic.co.uk)