
# AMBIC®

## AUTO FOAMER

### LV-REMOTE PRESSURE SWITCH - 2015

## OPERATING INSTRUCTIONS

MODE D'EMPLOI

ANLEITUNG

*leading best practice in livestock health management*

Fig. 1 – GENERAL LAYOUT OF AUTOFOAMER SYSTEM



Fig. 1B – TWIN PUMP AUTOFOAMER INSTALLATION – ENCLOSURE LAYOUT/ASSEMBLY

# Ambic Equipment Ltd - AutoFoamer Operating Notes

## Introduction

The AutoFoamer is designed to be installed using twin distribution tubes of 4mm bore (1/4" outside diameter) – the BLACK tube supplies compressed air; the GREY /BLUE tube supplies pressurised teat disinfectant chemical – both at a low pressure 3 - 5 psi (0.2 - 0.4 Bar) – see Fig 1. Chemical and air are supplied to the foaming applicators, which generate foam in the cup whilst the trigger is held depressed.

**WARNING** – Use only the GREY/BLUE tube for the CHEMICAL – the Black tube may be attacked by some chemicals.

## Safety

The AutoFoamer Foam Teat Dipping System is designed exclusively for use in milking installations. Any application outside the use described in this operating manual will be taken to be not in accordance with the intended purpose. The manufacturer/supplier will not be held responsible for any losses arising as a result of such use. The user will take full responsibility for use. USE IN ACCORDANCE WITH THE INTENDED PURPOSE ALSO INCLUDES COMPLYING WITH THE OPERATING MANUAL AND THE CONDITIONS FOR INSPECTION AND MAINTENANCE.

## Installation (refer to Fig.1 & 1B, or 1C)

**Electrical Installation** - should ideally be carried out by a Qualified Electrician – EXCEPT where a suitable waterproof plug outlet is available. THIS ELECTRICAL SUPPLY SHOULD BE PROTECTED BY AN APPROPRIATE EARTH LEAKAGE CIRCUIT BREAKER. IT IS RECOMMENDED THAT THE ELECTRICAL SUPPLY SHOULD BE SWITCHED OFF AUTOMATICALLY WHEN THE VACUUM PUMP IS NOT IN OPERATION.

**Pressurised Chemical** – is sucked up from a chemical container by peristaltic pump(s). The enclosure housing the pump(s) of the AutoFoamer is rated at IP40 (IEC 60529), the power unit (A) must therefore be positioned outside of the milking parlour in a dry location and close to a 220 –240V AC power socket (Power consumed: 1 pump = 120W ; 2 pumps = 240W). The 3-core power supply cord MUST remain accessible when the unit has been installed and MUST be connected to Earth.

The large enclosure (A1) should be positioned no more than 2 metres vertically above the floor and no more than 3 metres distant from the chemical container(s). For efficient and reliable operation the maximum vertical distance from pump to distribution line should be less than 1 metre. Use the drilling template (page 31) to mark out screw positions and fix the enclosure to a suitable flat wall, using the screws and plastic plugs supplied.

**NOTE:** Allow sufficient space to the Right of the enclosure for mounting the **Pressure Reservoir Bottle (W)**, which **MUST** be in an UPRIGHT position, using a single screw (see Fig. 1B or 1C).

**The small enclosure (A2) housing the Pressure Switch** (connected to the large enclosure by 4 metres of armoured cable) should be screwed to a flat vertical surface (e.g. Wall) some 1-4 metres distant from the large enclosure but as close as possible to the highest point of the distribution tubing. Use Cable clips supplied to secure cable, but DO NOT SHORTEN armoured cable.


**The chemical container(s) (D)** should be placed securely on the floor and the Inlet Tube(s) (C) from the pump enclosure connected to the Inlet nipple(s) of the peristaltic pump(s) – each PVC inlet tube should be warmed before pushing over the peristaltic pump nipple and securing with the stainless steel hose clips supplied. Make sure that clips are fitted so that the barb of the fitting on the pump is central in the clip and do NOT over tighten the clip (you risk breaking the plastic nipple). Drop the filter end of the inlet tube(s) into the chemical container(s).

**Low Pressure Air** is supplied, via the Regulator (B), from a DRY compressed air supply of 2-6 Bar (30-90 psi), using suitable connection fittings (not supplied). The enclosure containing the Regulator is sealed to IP55 and may be situated in any convenient, easily accessible location. The Non-Return Valve (V) is installed in the outlet line to protect the compressor in the event of any Non-Return Valve failure elsewhere in the system. During the initial system setting up process, it may be necessary to adjust the pressure regulator setting to obtain ideal operating conditions for efficient foam production.

continued . . .

Installation (continued)


**Fig. 1C – SINGLE PUMP AUTOFOAMER INSTALLATION – ENCLOSURE LAYOUT/ASSEMBLY**

**Twin Distribution tubes (Black & Grey /Blue)** - Routing of the distribution tubes (E) will depend on the parlour design and can be either at high level, or low level. In order to ensure efficient operation of the foaming applicators route the distribution line so it is as short as possible.

**Fig. 2 –TWIN TUBE FIXING BRACKETS – ASSEMBLY/FIXING**

Distribution tubes can be routed at high level, along the centre of the parlour, or at low level beneath the cow standing kerb. **Twin Tube Fixing Brackets (J)** are supplied and used to secure the main distribution tubing to suitable stallwork, using a single cable tie (**K**) around the tube to clamp the two halves of the brackets together (see Fig. 2). Alternatively, these Fixing Brackets may be affixed to a flat surface using a suitable screw (not supplied) through the central hole.

**T-pieces (F) and other distribution line fittings** are "push-in" and can be released by holding the flange in and pulling out the tubing – see diagrams fig. 4 – it is important that tubes are pushed fully into fittings as otherwise air (or liquid) will leak out from the tube/fitting.

**Foam Applicators** should be positioned such that all milking points can be reached comfortably without over-stretching the coils. Extension kits are available to enable as many Applicators as may be required to be installed.

**WARNING: Do NOT over-stretch the coils**, as this will place undue strain on the tubing and may result in the connectors, or manifold, leaking or becoming damaged. It is therefore NOT recommended to allow the applicator to hang from the coil when not in use, as this will result in over-stretching of the coil and possible damage to the applicator through contact with the floor.

Fig. 3 – ASSEMBLY OF AUTOFOAMER FOAM APPLICATORS



Fig. 4 – USING “PUSH-FIT” CONNECTIONS

## Ambic Equipment Ltd - AutoFoamer Operating Notes

### Installation (continued)

The **Plastic Hanging Hook** (see fig. 3, H) may be fitted to a suitable horizontal tube using 2 cable ties (K), or can be drilled with 2 holes and fitted (by 2 screws or bolts) to a flat surface. It is so designed that the cup of the applicator can rest in the Plastic Hanging Hook.

**Support Brackets** are supplied for each drop point; these serve both to hold the two T-pieces (F) and to support the coil ends to avoid a direct pull on the connections. The Support Brackets (L) are designed to be fitted to a horizontal tube using 2 cable ties (K) around the tubing (see fig. 3). If a suitable horizontal tube is not available, then Support Brackets should be fitted by alternative means.

The 2 T-pieces (F) of the distribution line locate into either side of the Support Bracket (L) and are secured by two cable ties (see fig. 3 K). The tubes from the twin tube coil (G) are pushed firmly into the respective T-pieces (F) and are twisted slightly to rest in the guides; the tubes are then clamped in place by fitting the final section of the bracket and securing it with a M4 bolt and lock-nut or a single cable tie (K). Although it is not essential that the Black and Blue tails of drop coils are connected to specific distribution lines, it is preferable that the BLACK tube is continuous across all drops and applicators.

**Non-Return Valves** (see fig. 1 V) - In this version of the AutoFoamer discreet Non-Return Valves are used. 2 Non-Return Valves are supplied with each Drop Point; these valves should be sited as indicated in fig. 5.


- 1 Non-Return Valve in the Air Supply tube
- 1 Non-Return Valve in the Liquid Supply tube

The valves should both be sited immediately before the Drop Coil Support Bracket with the GREY parts (arrows) pointing towards the Foam Applicator.

The Non-Return Valves should be pushed well in to fittings (typically 16mm minimum) to ensure proper sealing with the internal O-ring. Valves may be removed from the fittings by pushing in the collet then pulling out the tubing (as per other push-fit fittings used with this equipment).

Each Drop Point must have its own pair of Non-Return Valves – hence the use of the T-pieces in the distribution line where more than one Drop Point is fitted.

**WARNING –** the system will NOT produce foam correctly if the Non-Return Valves are missing or incorrectly fitted.

Fitting of the Foam Applicator to the Drop Coil is by way of the Manifold (M) which is delivered ready-fitted to one end of the Drop Coil. The Manifold is screwed into the rear of the trigger assembly of the Foaming Applicator and sealed with the O-ring which is fitted on the protruding end of the Manifold before screwing carefully into the rear end of the gun body. Use the AmbiSpanner supplied to tighten the manifold hand-tight to seal.

**WARNING - DO NOT ATTEMPT TO REMOVE THE MANIFOLD FROM THE COIL –** you risk causing permanent irreparable damage to the manifold (replacement drop coils are supplied with Manifold ready-fitted).

## Ambic Equipment Ltd - AutoFoamer Operating Notes

### Initial Setting Up

Leave the Foam Applicator at the far end of the distribution line disconnected and, in its place, fit a short length of tube resting in a suitable container in order to collect any chemical discharged.

Switch on the power at the electric socket and on the front of the pumping unit. Power Unit (**A1**) will be activated and liquid will start to be drawn up through the pump(s). Allow unit to run until liquid has passed through the pump(s) and is flowing out in a continuous stream from the end of the distribution pipe (this may take a few minutes on a system with a long distribution line). At this point, turn Off the unit and fit the final applicator. Briefly remove the plug from the T-piece below the small enclosure (**A2**) replacing it as soon as liquid runs out, to ensure liquid has reached the pressure switch.

Switch the pump unit back On and, when the Pressure Reservoir Bottle (**W**) has been filled to a level of about 30mm, and a liquid pressure of ~3psi (0.2 Bar) is reached, the pressure switch will automatically switch off the peristaltic pump(s), which will not operate until a gun is used and pressure drops. If pumps fail to stop within 1 minute (or run continuously) then check the distribution line system for leaks.

Switch on the air supply, adjust the pressure regulator initially to 3 – 4 psi and check the quality of foam by depressing the trigger on each applicator in turn. Initially the foam may not be produced until a continuous stream of liquid has reached the base of the foaming cup; priming each applicator should take no more than a few seconds and is accomplished by holding the trigger fully depressed until foam fills the applicator cup.

Foam Quality/Wetness may be adjusted by regulating the relative liquid to air pressure ratio. The unit will have been preset in the factory to a near optimum liquid pressure and it is therefore preferable to adjust the air pressure to arrive at the desired foam quality. Increasing Air pressure will make the foam more “dry” (very frothy and tending towards larger air bubbles), whilst decreasing air pressure will both reduce the rate of production of foam and make the foam “wetter” ultimately to the extent where foam collapses very quickly to leave fluid in applicator cup. The difference in air pressure between these extremes typically covers the pressure range of 2 – 6 psi (with liquid pressure set at 3 psi). Generally the Air pressure should be within  $\pm 2$  psi of Liquid pressure.

### OPERATION

When ready to dip, depress the trigger until foam is level or just above the rim of Applicator (F) and raise fully onto teat. Typically a full cup will dip two teats.

WARNING – Refill chemical container(s) BEFORE EMPTY and SWITCH OFF PUMP UNIT whilst refilling. In the event that chemical supply is exhausted, it will be necessary to re-prime system as described in the section “Initial Setting Up” above.

WHEN MILKING IS COMPLETED SWITCH OFF ELECTRIC POWER SUPPLY and rinse out cups with warm water.

**WARNING – AFTER MILKING IS COMPLETE ALWAYS ISOLATE MAINS POWER FROM PUMPS (this action is essential to avoid possible loss of chemical through leakage and damage to pump motors).**

### SPECIFICATIONS

Power Source – Liquid	– Electric Pumping unit	• 220-240Vac 50Hz 120/240W
	– Fuses	• Single Pump T0.4A • Twin Pump T1.0A
Power Source – Air	– Compressed Air	• 5 -10L/min at 3 - 7psi (0.2 - 0.5 Bar)
Maximum No. of Applicators per Pumping Unit		• 10
Maximum No. of Applicators operating simultaneously		• 2
Maximum Length of Distribution Tubing (Black or Grey/ Blue)		• 45 Metres (150 ft)
Typical Chemical Consumption per Applicator		• 1.5 – 2.5 mL/second
Operating Temperature		• 5 – 40 °C

## Ambic Equipment Ltd - AutoFoamer Operating Notes


### MAINTENANCE

**DAILY** - Applicator Cups should be rinsed out with clean warm water to remove any hairs, dirt, etc. immediately following completion of milking – this should avoid the build up of deposits and prevent chemical solidifying in the gauze. If foam production becomes slow and/or difficult then the Foaming Applicator cup may be easily dismantled (see fig. 6 below) for cleaning as follows:-

- Unscrew the cup assembly (N) and remove from the end of the lance (P) by gripping the outside of the cup and turning anti-clockwise.
- Carefully lift off and retain the large O-ring (S), then lift out the 3 Filter Gauzes (T) from the cup/retainer (R) at the end of the lance.

The gauzes may be cleaned by immersion in hot water and any stubborn material removed with a soft brush before re-assembling in the cup with O-ring on top.

**Fig. 6 – DISMANTLING FOAM APPLICATOR FOR CLEANING**

Use **ONLY** a soft dry cloth to clean the electrical enclosure housing when necessary – **NEVER** use a hose. If the supply cord is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified persons in order to avoid a hazard. Always disconnect the supply cord from the mains power supply **BEFORE** removing the enclosure lid.

### ANNUAL MAINTENANCE

In order to preserve the efficiency and reliability of the AutoFoamer system, we recommend that the following parts be replaced **EVERY 1 – 2 YEARS**, depending on number of cows dipped.

**Peristaltic pump tubes** – for **Single Pump** units use kit **AAF/022**.  
- for **Twin Pump** units use kit **AAF/023**.

**Pressure Switch** – use kit **AAF/013**.

## Ambic Equipment Ltd - AutoFoamer Operating Notes

### PROBLEMS & FAULT FINDING

**WARNING - Always disconnect the supply cord from the mains power supply BEFORE removing the enclosure lid.**

Problem	Likely Cause	Probable Solution
No foam produced when trigger squeezed (ALL units)	No Air Supply present  No liquid present	Switch on Air supply Check for leaking connections on distribution tubing system. Switch on pumping unit; check pump turns when gun trigger depressed. Check/replenish chemical and prime system, if necessary.
No foam produced when trigger squeezed (at one unit only)	Air supply and/or liquid not reaching foaming applicator.	Leaking connections, or blocked tubing – check and remedy. Trigger valve stuck or blocked – unscrew valve and check, or replace with AAF/015. Foaming cup dirty or blocked – clean to remove dirt from gauzes, or replace with new cup.
Foam very dry, or a few large bubbles only produced when trigger squeezed.	No liquid supply or restriction in liquid distribution pipe work.	Check chemical container not empty. Check inlet tubing is always full of liquid and no leaks at connection to pump(s) Check distribution tubing (Grey/Blue) for continuity of liquid supply. Check Pump(s) running when trigger squeezed and held open. Check Air pressure not set too high at regulator – reduce pressure to 3-4 psi, if necessary.
Foam very wet, or only liquid rises in cup when trigger squeezed.	No air supply, or air distribution pipe work restricted/blocked.  Foaming cup air feed blocked/restricted.	Check Air supply switched On. Check for leaking connections on air distribution tubing (Black) system. Check Air pressure not set too Low at regulator. Unscrew applicator cup from lance; remove gauze disks; turn Off Pumping Unit; depress applicator trigger and check for Air flowing from hole in centre of Cup. Check that Air is reaching Applicator at gun body entry. Clean/ unblock cup and tubing, as appropriate.
No liquid reaching applicator drops.	Pump(s) not turning when triggers held squeezed.  Pump(s) turning but not pumping liquid.	Check fuses and replace if necessary. Check for blockages in Inlet and Outlet pipes. Possible Pressure Switch failure (inside small enclosure – <b>A2</b> ) - consult engineer. Check for leakage (or blockage) in Inlet pipe.  Possible failure of peristaltic pump tubing – consult engineer for replacement.

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

### AutoFoamer – LIQUID PRESSURE SETTING

During the initial installation and if it becomes necessary to replace the Pressure Switch, it may be necessary to adjust the liquid pressure slightly.

The Pressure Switch in all units despatched after 1 January 2015 is controlled via a transformer at 12 VAC. This switch allows the liquid pressure level to be adjusted in the range 2-5 psi (0.2-0.3 Bar) by means of a setscrew on the switch.

1. Turn off the mains power at the switch and disconnect the unit from the Mains Power supply.
2. Unscrew the 4 screws and lift the lid of the small enclosure (**A2**) to allow access to the inside of the enclosure; the pressure switch is arrowed in the photograph opposite.
3. Use a 5/64" (2mm) Hexagon Key to turn the recessed screw on the switch at the opposite end to the electrical connections (see next photo below opposite).
4. When viewed from above, turning the screw **CLOCKWISE** will **INCREASE** PRESSURE; **ANTICLOCKWISE** will **DECREASE** liquid Pressure. (see diagram below).

5. Each full turn of the adjusting screw increases/decreases the pressure by ~1 psi (~0.1 Bar).

It is recommended that **NO MORE THAN 3 TURNS** of the adjusting screw be made in either direction. **IF THE SCREW BECOMES LOOSE TO TURN, DO NOT TURN** any further as you risk causing permanent damage to the switch.

6. Remove the Hexagon Key and close the lid before re-connecting the unit to Mains Power to check if further adjustments are necessary.

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**Ambic Equipment Limited,