

Product Technical Data Sheet

EasiFoamer™Applicator & Extension Kit AAF/004 & AAF/014





Properties

Shield, Lance and Hook: Polypropylene

Gun Body & Lever and Connectors: Polypropylene

Mesh: Stainless Steel

Retractable Coil: Nylon blend 'o' rings: Silicone and Viton

Cable Ties: Nylon

Dimensions: Applicator = $18 \times 50 \text{ cm}$ bag Extension Kit = $18 \times 17 \times 33 \text{ cm}$

Net weight: Applicator = 142 g Extension Kit = 759 g

Background

Pre- and post-milking dairy hygiene is recommended for the pro-active control of Mastitis-associated pathogens. Whilst post-milking treatments have been effective at reducing the level of contagious mastitis, pre-milking treatment of cows' teats has been found to reduce the risk of contamination and new infections, particularly those caused by environmental pathogens. Foam is particularly effective at loosening and lifting dirt and faecal matter from the teat as it is better at penetrating to the skin surface than liquid. The use of foaming chemicals also reduces chemical consumption when compared to using liquid dips or sprays.



Product Technical Data Sheet

EasiFoamer™ applicators use special patented foaming mesh inserts to generate optimum foam structure ensuring that the appropriate amount of chemical is applied to the teat each time. They also facilitate the correct bubble structure that is essential for agitation of the teat surface in order to lift dirt and debris.

Features

- Especially designed for fast foam application.
- Minimises labour compared to manual filling of foaming dip cups.
- Easy to dismantle for rapid cleaning. Handy bracket enables applicator to be hung up, reducing wear and tear on the applicator and coil.
- Extension kits enable expansion of the number of drop coils on the EasiFoamer™ system.
- Patented system of foaming meshes to generate foam of an optimum quality (Patent EP2124524).
- Applies just the right amount of chemical to the teat requiring a single towel wipe to remove foaming chemicals.
- The following parts (in addition to those shown above) can be ordered separately as product codes:
 - AAF/004 AutoFoamer™ Lance Assembly.
 - AAF/006 AutoFoamer™ Coil with Manifold.
 - AAF/007 AutoFoamer™ Filter Disk Pack.
 - o AAF/008 AutoFoamer™ 'O'-Ring Pack (3 large & 3 small).
 - AAF/015 AutoFoamer™ Trigger Valve Kit.
 - AAF/016 AutoFoamer™ Manifold Repair Kit.
 - AAF/017 AutoFoamer™ Applicator Hanging Hook.
 - AAF/018 AutoFoamer™ Dropper Support Bracket.
 - AAF/019 AutoFoamer™ T Piece (Twin Pack).
 - ATS/405 Trigger Valve Service Kit.
 - ATS/406 Retractable Coil.
 - ATS/407 T Connector with Nuts.
 - o ATS/408 Straight Connectors with Nuts (pack of 2).
 - ATS/409 Connector Nuts (pack of 6).
 - ATS/410 6m Flexible Delivery line tubing (included in complete extension kits).
 - ATS/411 Cable Ties (Pack of 10).
 - ATS/438 Repair Kit for Ambic Drop Coil.
 - ATS/461 Straight Connector.
 - AJS/2055 AmbiSpanner™.

Country of Origin: GB EU Tariff Code: 84349000

REACH compliance is not applicable.