

AMBIC®

...designed for today's Dairy Farmer

Choose the right dip cup for your application

Model/ Feature	Standard Dipper	No.1 Non Return	Mini Dipper	Twin Dipper	Uni Dipper	Multi Foamer
Non Return		✓	✓	✓	✓	
Standard Dips	✓	✓	✓		✓	
Barrier Dips				✓	✓	
Foaming Dips						✓
15ml Cup size			✓			
30ml Cup size	✓	✓		✓	✓	✓
Cows	✓	✓	✓	✓	✓	✓
Goats					✓	
Sheep					✓	
Small teated animals			✓			✓
DipMizer compatible	✓	✓	✓	✓	✓	

Part numbers:

ADM/101-B	DipMizer - single pack - Blue
ADM/103-B	DipMizer - pack of three - Blue
ADM/030-5	DipMizer brush only - pack of five

YOUR LOCAL AMBIC DEALER:



Ambic Equipment Limited
 Avenue Four, Station Lane,
 Witney, Oxfordshire,
 OX28 4XT. England
 Tel: +44 (0) 1993 701936/7
 Fax: +44 (0) 1993 779039
 Email: sales@ambic.co.uk



SHEET No. DM6402 25.06.08

© Ambic Equipment Limited 25.06.08

AMBIC®

DipMizer™

From the world's leading dip cup supplier



- Removes excess teat dip
- Easily removed for regular routine cleaning
- Easy cleaning means less risk of contamination
- Retrofits all Ambic (non foaming) dip cups
- Low density bristles remove just enough chemical, but not too much!

Visit the Ambic web site: www.ambic.co.uk

leading best practice in livestock health management

AMBIC®

How the DipMizer works

- The traditional dip cup applies teat dip by immersing the teat in chemical and coating it with teat dip.
- Adding the DipMizer™ device adds a 'sweeping' step to the dipping routine in order to remove excess chemical from the cow's teat as it is withdrawn from the dip cup – saving you money on unnecessary teat dip consumption.
- Removing too much teat dip can be counter-productive in that udder health could be compromised by leaving too little coating on the teat. The DipMizer™ is designed with special bristles that remove just enough chemical, but not too much so that the teat is well protected.
- Always remove the DipMizer™ brush attachment at the end of milking and sanitise it in a cleaning solution to reduce the risk of contamination.
- Use a fresh clean brush attachment for each milking as recommended in an evaluation study undertaken by the Institute for Animal Health, Compton, UK.

AMBIC®

Use DipMizer with these DipCups

• **Original No1. Non Return dip cup™**

The world's first non return dip cup that leads the market in quality and remains the world's no.1 selling dip cup. 30ml dip cup capacity with non return design which prevents contaminated chemical returning to the dip cup bottle.

• **MiniDipper™**

An alternative to the popular Original No. 1 Non Return dip cup with a 15ml dip cup capacity to suit smaller teated cows.

• **TwinDipper™**

An alternative version to the Original No.1 Non- Return dip cup designed for ease of application of thicker, barrier type or 'film forming' dips which have higher viscosity. Twin application tubes reduce the effort required to dispense the dip.

• **UniDipper™**

'Low profile' design suitable for dipping cow, goat and sheep teats or for farmers who prefer a 'side action' dipping routine, such as in a stanchion barn. Unique DipWell™ feature to prevent spillage and choice of single or twin tube application to suit most dips.

• **StandardDipper™**

The world's very first dip cup with large angled cup for ease of application and maximum teat coverage.

