

## Pump Head Diaphragm

<p style="text-align: center;"><b>ATS/456</b></p> 	<p style="text-align: center;"><b>ATS/456-V</b></p> 
<p>Diaphragm: EPDM Stud: Bonded Mild Steel Dimensions: 8 x 8 x 3.5 cm Net Weight: 29 g</p>	<p>Diaphragm: Viton® Stud: Bonded Mild Steel Dimensions: 8 x 8 x 3.5 cm Net Weight: 42 g</p>
<p>Used in: JetStream™ Classic EasiDipper PeraSpray™</p>	<p>Used in: Viton JetStream™ Viton EasiDipper OptiFlow</p>
<p>Suitable for use with*: Iodine Chlorhexidine Peracetic Acid</p>	<p>Suitable for use with*: Iodine Chlorhexidine Lactic Acid Chlorine Dioxide</p>
<p><b>Not</b> suitable for use with*: Lactic Acid Chlorine Dioxide Sodium Hypochlorite</p>	<p><b>Not</b> suitable for use with*: Peracetic Acid Sodium Hypochlorite</p>

**\*Please note that the chemical listings are not comprehensive.**

**Each diaphragm pump contains 2 diaphragms. We recommend that both should be replaced at the same time to avoid excessive wear and tear.**

Manufactured in the UK  
HS Commodity Code: 84349000